

Release Notes: MATLAB Plug-in for Teledyne Lumenera camera model's support

Requirements:

- The Microsoft Visual Studio 2015 redistributable files are required and be downloaded directly from the Microsoft Website.
- The Teledyne [LuCam Software and Software Development Kit](#) is required for both camera API and camera device drivers. (Minimum requirement is v6.8.2)
- This plug-in has been tested with MATLAB R2016a/b, R2017a/b, R2018a/b, R2019a.

Changes from 2.1.0 to 2.1.1

Command line interface

- Adding command `LucamGetSnapshot()` and `LucamSetSnapshot()` to interface to Teledyne Lumenera API structure.

Imaqtools interface

- Add support for R2018a/b, R2019a.
- Add AE (Auto Exposure) property. This allow to enable or disable camera Auto Exposure property in video mode.

Changes from 1.X.X to 2.1.0

Command line interface

- Adds the ability to capture video or snapshot raw frame and convert it later.
- Raw frame is a structure containing information about frame format, conversion parameters and image data.
- `LucamTurboxxxxxx` functions added to capture raw frames into memory and converts it later when you return frames into Matlab.
- `LucamGetProperty` and `LucamGetPropertyRange` will return a data structure containing all the information.
- `LucamSetProperty` will use the data structure output from `LucamGetProperty`.
- The `.m` files for all commands been updated to throw `MException` on error.
- Overall code re-factoring.
- Started the process of altering the implementation mechanism of commands for easier support maintenance.

Imaqtools interface

- Added support of USB 3.0 cameras.
- Added more frame formats selections.
- Changed all Properties names to be similar to the API.
- Properties will be populated based on what is supported by the camera model.
- Added a property read-back function.